

Product Datasheet

Human ACTg2 ELISA Kit (Ready to Use) (orb549688)

Description	Human Actin Gamma 2, Smooth Muscle (ACTg2) Sandwich ELISA KitReady to Use: Pre-diluted detection rea...
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates or other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	P63267
Sensitivity	0.053ng/mL
Entrez	72
Expiration Date	Please enquire.