

Product Datasheet

Human ARPC2 ELISA Kit (Ready to Use) (orb549689)

Description	Human Actin Related Protein 2/3 Complex Subunit 2 (ARPC2) Sandwich ELISA KitReady to Use: Pre-dilute...
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates or other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	O15144
Sensitivity	0.068ng/mL
Entrez	10109
Expiration Date	Please enquire.