

Product Datasheet

Human ATF6 ELISA Kit (orb549698)

Description	Human Activating Transcription Factor 6 (ATF6) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Types	tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Sensitivity	0.049ng/mL
Expiration Date	Please enquire.