

Product Datasheet

Human ADORA2A ELISA Kit (Ready to Use) (orb549710)

Description	Human Adenosine A2a Receptor (ADORA2A) Sandwich ELISA KitReady to Use: Pre-diluted detection reagent...
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates or other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	P29274
Sensitivity	0.060ng/mL
Entrez	135
Expiration Date	Please enquire.