

Product Datasheet

Human beta 1 Adrenergic Receptor ELISA Kit (Ready to Use) (orb549730)

Description	Human Adrenergic Receptor Beta 1 (beta 1 Adrenergic Receptor) Sandwich ELISA KitReady to Use: Pre-di...
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Sensitivity	0.058ng/mL
Expiration Date	Please enquire.