

Product Datasheet

Human ADM2 ELISA Kit (Ready to Use) (orb549736)

Description	Human Adrenomedullin 2 (ADM2) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for con...
Reactivity	Human
Range	78.1-5000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	Q7Z4H4
Sensitivity	26.2pg/mL
Entrez	79924
Expiration Date	Please enquire.