

Product Datasheet

Human ALDH7A1 ELISA Kit (Ready to Use) (orb549756)

Description	Human Aldehyde Dehydrogenase 7 Family, Member A1 (ALDH7A1) Sandwich ELISA KitReady to Use: Pre- dilut...
Reactivity	Human
Range	0.781-50ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Sensitivity	0.29ng/mL
Expiration Date	Please enquire.