

Product Datasheet

Human ALDOB ELISA Kit (Ready to Use) (orb549760)

Description	Human Aldolase B, Fructose Bisphosphate (ALDOB) Sandwich ELISA KitReady to Use: Pre-diluted detectio...
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P05062
Sensitivity	0.063ng/mL
Entrez	229
Expiration Date	Please enquire.