

Product Datasheet

Human TTPa ELISA Kit (Ready to Use) (orb549770)

Description	Human Alpha Tocopherol Transfer Protein (TTPa) Sandwich ELISA KitReady to Use: Pre-diluted detection...
Reactivity	Human
Range	62.5-4000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates or other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Sensitivity	25.3pg/mL
Expiration Date	Please enquire.