

Product Datasheet

Human a1AGP ELISA Kit (Ready to Use) (orb549771)

Description	Human Alpha-1-Acid Glycoprotein (a1AGP) Sandwich ELISA KitReady to Use: Pre-diluted detection reagen...
Reactivity	Human
Range	3.12-200ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	P02763
Sensitivity	1.27ng/mL
Entrez	5004
Expiration Date	Please enquire.