

Product Datasheet

Human a1ACT ELISA Kit (Ready to Use) (orb549772)

Description	Human Alpha-1-Antichymotrypsin (a1ACT) Sandwich ELISA KitReady to Use: Pre-diluted detection reagent...
Reactivity	Human
Range	3.12-200ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P01011
Sensitivity	1.42ng/mL
Entrez	12
Expiration Date	Please enquire.