

Product Datasheet

Human a2ML1 ELISA Kit (Ready to Use) (orb549779)

Description	Human Alpha-2-Microglobulin Like Protein 1 (a2ML1) Sandwich ELISA KitReady to Use: Pre-diluted detec...
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	A8K2U0
Sensitivity	0.067ng/mL
Entrez	144568
Expiration Date	Please enquire.