

Product Datasheet

Human aMACR ELISA Kit (Ready to Use) (orb549787)

Description	Human Alpha-Methylacyl Coenzyme A Racemase (aMACR) Sandwich ELISA KitReady to Use: Pre-diluted detec...
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates, cell lysates or other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q9UHK6
Sensitivity	0.059ng/mL
Entrez	23600
Expiration Date	Please enquire.