

Product Datasheet

Human ANGPTL1 ELISA Kit (Ready to Use) (orb549821)

Description	Human Angiopoietin Like Protein 1 (ANGPTL1) Sandwich ELISA KitReady to Use: Pre-diluted detection re...
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	O95841
Sensitivity	0.055ng/mL
Entrez	9068
Expiration Date	Please enquire.