

Product Datasheet

Human ACE2 ELISA Kit (Ready to Use) (orb549832)

Description	Human Angiotensin I Converting Enzyme 2 (ACE2) Sandwich ELISA KitReady to Use: Pre-diluted detection...
Reactivity	Human
Range	15.6-1000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q9BYF1
Sensitivity	7.7pg/mL
Entrez	59272
Expiration Date	Please enquire.