

Product Datasheet

Human ANXA11 ELISA Kit (Ready to Use) (orb549842)

Description	Human Annexin A11 (ANXA11) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for conven...
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates or other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	P50995
Sensitivity	0.061ng/mL
Entrez	311
Expiration Date	Please enquire.