

Product Datasheet

Human AQP9 ELISA Kit (Ready to Use) (orb549888)

Description	Human Aquaporin 9 (AQP9) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for convenie...
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	O43315
Sensitivity	0.054ng/mL
Entrez	366
Expiration Date	Please enquire.