

Product Datasheet

Human ALOXE3 ELISA Kit (Ready to Use) (orb549889)

Description	Human Arachidonate Lipoxygenase 3 (ALOXE3) Sandwich ELISA KitReady to Use: Pre-diluted detection rea...
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates, cell lysates or other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q9BYJ1
Sensitivity	0.114ng/mL
Entrez	59344
Expiration Date	Please enquire.