

Product Datasheet

Human ABCC9 ELISA Kit (Ready to Use) (orb549938)

Description	Human ATP Binding Cassette Transporter C9 (ABCC9) Sandwich ELISA KitReady to Use: Pre-diluted detect...
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates, cell lysates or other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	O60706
Sensitivity	0.132ng/mL
Entrez	10060
Expiration Date	Please enquire.