

Product Datasheet

Human ATP2A2 ELISA Kit (Ready to Use) (orb549941)

Description	Human ATPase, Ca++ Transporting, Cardiac Muscle, Slow Twitch 2 (ATP2A2) Sandwich ELISA KitReady to U...
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates, cell lysates or other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	P16615
Sensitivity	0.126ng/mL
Entrez	488
Expiration Date	Please enquire.