

Product Datasheet

Human ATP2B2 ELISA Kit (Ready to Use) (orb549943)

Description	Human ATPase, Ca++ Transporting, Plasma Membrane 2 (ATP2B2) Sandwich ELISA KitReady to Use: Pre- dilu...
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates, cell lysates or other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Sensitivity	0.112ng/mL
Expiration Date	Please enquire.