

Product Datasheet

Human ATP6AP2 ELISA Kit (Ready to Use) (orb549946)

Description	Human ATPase, H ⁺ Transporting, Lysosomal Accessory Protein 2 (ATP6AP2) Sandwich ELISA Kit Ready to Use...
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability tissue homogenates or other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Sensitivity	0.112ng/mL
Expiration Date	Please enquire.