

Product Datasheet

Human ATP4b ELISA Kit (Ready to Use) (orb549949)

Description	Human ATPase, H ⁺ /K ⁺ Exchanging Beta Polypeptide (ATP4b) Sandwich ELISA KitReady to Use: Pre-diluted ...
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates or other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Sensitivity	0.059ng/mL
Expiration Date	Please enquire.