

Product Datasheet

Human ATP1a1 ELISA Kit (orb549950)

Description	Human ATPase, Na ⁺ /K ⁺ Transporting Alpha 1 Polypeptide (ATP1a1) Sandwich ELISA Kit
Reactivity	Human
Range	78.125-5000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates or other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Sensitivity	26.6pg/mL
Expiration Date	Please enquire.