

Product Datasheet

Human ATP1b4 ELISA Kit (Ready to Use) (orb549957)

Description	Human ATPase, Na ⁺ /K ⁺ Transporting Beta 4 Polypeptide (ATP1b4) Sandwich ELISA KitReady to Use: Pre-di...
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates, cell lysates or other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q9UN42
Sensitivity	0.106ng/mL
Entrez	23439
Expiration Date	Please enquire.