

Product Datasheet

Human ATG3 ELISA Kit (Ready to Use) (orb549964)

Description	Human Autophagy Related Protein 3 (ATG3) Sandwich ELISA KitReady to Use: Pre-diluted detection reage...
Reactivity	Human
Range	15.6-1000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates, cell lysates or other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Sensitivity	6.9pg/mL
Expiration Date	Please enquire.