

Product Datasheet

Human BIRC6 ELISA Kit (Ready to Use) (orb549977)

Description	Human Baculoviral IAP Repeat Containing Protein 6 (BIRC6) Sandwich ELISA KitReady to Use: Pre-dilute...
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates, cell lysates or other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q9NR09
Sensitivity	0.107ng/mL
Entrez	57448
Expiration Date	Please enquire.