


Product Datasheet

Human CEA / Carcino Embryonic Antigen ELISA Kit (orb54999)

Description	Human CEA / Carcino Embryonic Antigen ELISA Kit
Reactivity	Human
Range	312 pg/ml - 20,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: K35-A685
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	76795 MW
Uniprot ID	P06731
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum and plasma is around 400 pg/ml.
Expiration Date	Please enquire.



Human
CEA
PicoKine
ELISA Kit
standard
cu...