

Product Datasheet

Human bACE1 ELISA Kit (Ready to Use) (orb549995)

Description	Human Beta-Site APP Cleaving Enzyme 1 (bACE1) Sandwich ELISA KitReady to Use: Pre-diluted detection ...
Reactivity	Human
Range	0.25-16ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, cerebrospinal fluid or other biological fluids.
Sample Types	serum, plasma, cerebrospinal fluid and other biological fluids.
Uniprot ID	P56817
Sensitivity	0.11ng/mL
Entrez	23621
Expiration Date	Please enquire.