


Product Datasheet

Human Angiostatin Kringle 1-3/PLG ELISA Kit (orb55000)

Description	Human Angiostatin Kringle 1-3/PLG ELISA Kit
Reactivity	Human
Range	156 pg/ml - 10,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: V98-S355 (K1-3)
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA)
MW	90569 MW
Uniprot ID	P00747
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:1000, concentration in serum and plasma is around 1000 ng/ml.
Expiration Date	Please enquire.



Human
Angiostatin
K1-3
PicoKine
ELISA Ki...