

Product Datasheet

Human BLVRA ELISA Kit (Ready to Use) (orb550001)

Description	Human Biliverdin Reductase A (BLVRA) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents ...
Reactivity	Human
Range	0.625-40ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P53004
Sensitivity	0.282ng/mL
Entrez	644
Expiration Date	Please enquire.