

Product Datasheet

Human BMPR1A ELISA Kit (Ready to Use) (orb550016)

Description	Human Bone Morphogenetic Protein Receptor 1A (BMPR1A) Sandwich ELISA KitReady to Use: Pre-diluted de...
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates, cell lysates or other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	P36894
Sensitivity	0.11ng/mL
Entrez	657
Expiration Date	Please enquire.