

Product Datasheet

Human BASP1 ELISA Kit (Ready to Use) (orb550020)

Description	Human Brain Abundant, Membrane Attached Signal Protein 1 (BASP1) Sandwich ELISA KitReady to Use: Pre...
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates, cell lysates or other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	P80723
Sensitivity	0.112ng/mL
Entrez	10409
Expiration Date	Please enquire.