

Product Datasheet

Human BAI1 ELISA Kit (Ready to Use) (orb550024)

Description	Human Brain Specific Angiogenesis Inhibitor 1 (BAI1) Sandwich ELISA KitReady to Use: Pre-diluted det...
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	014514
Sensitivity	0.065ng/mL
Entrez	575
Expiration Date	Please enquire.