

Product Datasheet

Human BCAR1 ELISA Kit (Ready to Use) (orb550026)

Description	Human Breast Cancer Anti-Estrogen Resistance 1 (BCAR1) Sandwich ELISA KitReady to Use: Pre-diluted d...
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P56945
Sensitivity	0.053ng/mL
Entrez	9564
Expiration Date	Please enquire.