

Product Datasheet

Human C4BPb ELISA Kit (Ready to Use) (orb550051)

Description	Human C4 Binding Protein Beta (C4BPb) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents...
Reactivity	Human
Range	0.468-30ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P20851
Sensitivity	0.17ng/mL
Entrez	725
Expiration Date	Please enquire.