

Product Datasheet

Human CALB ELISA Kit (Ready to Use) (orb550061)

Description	Human Calbindin (CALB) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for convenient...
Reactivity	Human
Range	12.5-800pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P05937
Sensitivity	5.7pg/mL
Entrez	793
Expiration Date	Please enquire.