

Product Datasheet

Human CAMKK2 ELISA Kit (Ready to Use) (orb550070)

Description	Human Calcium/Calmodulin Dependent Protein Kinase Kinase 2 (CAMKK2) Sandwich ELISA KitReady to Use: ...
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates or other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	Q96RR4
Sensitivity	0.114ng/mL
Entrez	10645
Expiration Date	Please enquire.