

Product Datasheet

Human CIP2A ELISA Kit (Ready to Use) (orb550091)

Description	Human Cancerous Inhibitor Of PP2A (CIP2A) Sandwich ELISA KitReady to Use: Pre-diluted detection reagent...
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q8TCG1
Sensitivity	0.057ng/mL
Entrez	57650
Expiration Date	Please enquire.