

Product Datasheet

Human CSN1 ELISA Kit (Ready to Use) (orb550143)

Description	Human Casein Alpha (CSN1) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for conveni...
Reactivity	Human
Range	1.56-100ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates, cell lysates, cell culture supernates, breast milk or other biological fluids.
Sample Types	tissue homogenates, cell lysates, breast milk, cell culture supernatants and other biological fluids.
Uniprot ID	P47710
Sensitivity	0.55ng/mL
Entrez	1446
Expiration Date	Please enquire.