

Product Datasheet

Human CDC25 ELISA Kit (Ready to Use) (orb550187)

Description	Human Cell Division Cycle Protein 25 (CDC25) Sandwich ELISA KitReady to Use: Pre-diluted detection r...
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P30304
Sensitivity	0.062ng/mL
Entrez	993
Expiration Date	Please enquire.