

Product Datasheet

Human CRABP2 ELISA Kit (Ready to Use) (orb550190)

Description	Human Cellular Retinoic Acid Binding Protein 2 (CRABP2) Sandwich ELISA KitReady to Use: Pre-diluted ...
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates, cell lysates or other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	P29373
Sensitivity	0.114ng/mL
Entrez	1382
Expiration Date	Please enquire.