

Product Datasheet

Human CEMP1 ELISA Kit (Ready to Use) (orb550191)

Description	Human Cementum Protein 1 (CEMP1) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for ...
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q6PRD7
Sensitivity	0.124ng/mL
Entrez	752014
Expiration Date	Please enquire.