

Product Datasheet

Human CBLN1 ELISA Kit (Ready to Use) (orb550195)

Description	Human Cerebellin 1 (CBLN1) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for conven...
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P23435
Sensitivity	0.055ng/mL
Entrez	869
Expiration Date	Please enquire.