

Product Datasheet

Human CCR5 ELISA Kit (Ready to Use) (orb550209)

Description	Human Chemokine C-C-Motif Receptor 5 (CCR5) Sandwich ELISA KitReady to Use: Pre-diluted detection re...
Reactivity	Human
Range	15.6-1000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates, cell lysates or other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	P51681
Sensitivity	6.5pg/mL
Entrez	1234
Expiration Date	Please enquire.