

Product Datasheet

Human XCL2 ELISA Kit (Ready to Use) (orb550212)

Description	Human Chemokine C-Motif Ligand 2 (XCL2) Sandwich ELISA KitReady to Use: Pre-diluted detection reagen...
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q9UBD3
Sensitivity	0.058ng/mL
Entrez	6846
Expiration Date	Please enquire.