

Product Datasheet

Human CHRNa1 ELISA Kit (Ready to Use) (orb550233)

Description	Human Cholinergic Receptor, Nicotinic, Alpha 1 (CHRNa1) Sandwich ELISA KitReady to Use: Pre-diluted ...
Reactivity	Human
Range	0.781-50ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates or other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	P02708
Sensitivity	0.28ng/mL
Entrez	1134
Expiration Date	Please enquire.