

Product Datasheet

Human CHRNa3 ELISA Kit (Ready to Use) (orb550234)

Description	Human Cholinergic Receptor, Nicotinic, Alpha 3 (CHRNa3) Sandwich ELISA KitReady to Use: Pre-diluted ...
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates or other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	P32297
Sensitivity	0.107ng/mL
Entrez	1136
Expiration Date	Please enquire.