

Product Datasheet

Human CHRNb2 ELISA Kit (Ready to Use) (orb550237)

Description	Human Cholinergic Receptor, Nicotinic, Beta 2 (CHRNb2) Sandwich ELISA KitReady to Use: Pre-diluted d...
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates or other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	P17787
Sensitivity	0.054ng/mL
Entrez	1141
Expiration Date	Please enquire.