

Product Datasheet

Human CCDC3 ELISA Kit (Ready to Use) (orb550295)

Description	Human Coiled Coil Domain Containing Protein 3 (CCDC3) Sandwich ELISA KitReady to Use: Pre-diluted de...
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	Q9BQI4
Sensitivity	0.126ng/mL
Entrez	83643
Expiration Date	Please enquire.